



College of Pharmacy-University of Mustansiriyah
5Th. Year- Practical advance pharmaceutical analysis – 2019 -

Experiment-3

Titration of the ascorbic acid (vitamin C) in tablets By pH meter used first and 2nd derivatives


Outcomes:

After completing this experiment, the student should be able to:

- Calibration of pH meter.
- Calculation of first and 2nd derivatives.
- Calculate the K_a of weak acid.



Introduction:

In this experiment students analyze a tablet of Vitamin C (or other drug containing this compound) and calculate the percentage of ascorbic acid, $H_2C_6H_6O_6$, present in it.



Vitamin C (ascorbic acid)

Vitamin C tablets contain ascorbic acid as the active ingredient; however, it is mixed with such fillers as starch which however does not obscure the endpoint. Ascorbic acid has two steps of hydrolysis (hydrogen ion), ($pK_{a1}=4.10$, $pK_{a2}=11.8$). Since the second one comes off only at a pH of 10 or 11 the reaction of ascorbic acid with sodium hydroxide will produce the acid salt, sodium hydrogen ascorbate, and not the normal salt:-



Titration curve of 0.1 M ascorbic acid using 0.1 M NaOH as titrant

When titrated using an indicator like Phenolphthalein or Bromothymol Blue; The NaOH solution will be **standardized** to determine its strength by reacting it with a very pure (**primary standard**) sample of an acid.



College of Pharmacy-University of Mustansiriyah
5Th. Year- Practical advance pharmaceutical analysis – 2019 -

Materials and Equipment:- pH meter, buratte, pipette, beker, Vitamin C tablet, Sodium hydroxide, Indicator.

Procedure:

Place a tablet of Vitamin C in an Erlenmeyer flask, add about 50 mL of warm distilled water, and crush the tablet with a glass stirring rod. Cool down the solution. Add two or three drops of indicator to the flask and titrate with sodium hydroxide solution until the endpoint is reached. Make two independent titrations at least.

Result calculation:

As seen in the titration curve, vitamin C behaves as mono-protic acid when titrated using Phenolphthalein as indicator. This means that the number of moles of NaOH in titration is equal to the number of moles of the ascorbic acid.

Calculate the mass of acid for two titrations and, finally, calculate the arithmetic average of these two results. Compare the result with the factory value.

Note: that since this is a quantitative lab. exercise all measurements to be used in calculations must be recorded to the proper number of significant digits.

Molecular mass of ascorbic acid is 176.13 g/mol.

Record all the numbers obtained and the name of indicator used, as well as the calculations made.

Procedure (pH-metric titration):

1. Place a tablet of Vitamin C in an Erlenmeyer flask, add about 50 mL of warm distilled water, and crush the tablet with a glass stirring rod. Cool down the solution.
2. Dilute the sample in your volumetric flask to the total volume of 100-150 mL with distilled water.
3. Place the beaker on magnetic stirrer insert the magnet in it, as well as pH-electrode. Ask the assistant to control the correctness of the installation and some advices.
4. Measure and note pH; Repeat this adding small portions of the titrant base (few drops each), noting also the actual volume of titrant added (total). Finish titration when pH exceeds 12.

Processing the results

Using a computer data-sheet, plot the titration curve obtained and determine as exactly as possible the volume of titrant corresponding to the inflection point.

For more aspiring students:

The preciseness of determination of the inflection point can be very improved if calculating the second derivative of the titration curve.



College of Pharmacy-University of Mustansiriyah
5Th. Year- Practical advance pharmaceutical analysis – 2019 -

Record all the numbers obtained, as well as the calculations made.
Note your observation concerning the comparison of your result with the factory value.

Excel Shortcut Keystrokes for the PC*	
<i>*Macintosh equivalents, if different, appear in square brackets</i>	
<u>TO ACCOMPLISH THIS TASK</u>	<u>TYPE THESE KEYSTROKES</u>
Alternate between displaying cell values and displaying cell formulas	Ctrl+` [⌘+`]
Calculate all sheets in all open workbooks	F9
Calculate the active worksheet	Shift+F9
Cancel an entry in a cell or formula bar	Esc
Complete a cell entry and move down in the selection	Enter [Return]
Complete a cell entry and move to the left in the selection	Shift+Tab
Complete a cell entry and move to the right in the selection	Tab
Complete a cell entry and move up in the selection	Shift+Enter
Copy a formula from the cell above the active cell into the cell or the formula bar	Ctrl+' (Apostrophe) [⌘+']
Copy a selection	Ctrl+C [⌘+C]
Copy the value from the cell above the active cell into the cell or the formula bar	Ctrl+Shift+'' (Quotation Mark) [⌘+Shift+'']
Cut a selection	Ctrl+X [⌘+X]
Define a name	Ctrl+F3 [⌘+F3]
Delete the character to the left of the insertion point, or delete the selection	Backspace [Delete]
Delete the character to the right of the insertion point, or delete the selection	Delete [Del]
Displays the Insert Function dialog box	Shift+F3
Displays Key Tips for ribbon shortcuts	ALT
Edit a cell comment	Shift+F2
Edit the active cell	F2 [None]
Edit the active cell and then clear it, or delete the preceding character in the active cell as you edit the cell contents	Backspace [Delete]
Enter a formula as an array formula	Ctrl+Shift+Enter
Fill down	Ctrl+D [⌘+D]
Fill the selected cell range with the current entry	Ctrl+Enter [None]
Fill to the right	Ctrl+R [⌘+R]
Format cells dialog box	Ctrl+1 [⌘+1]
Insert the AutoSum formula	Alt+= (Equal Sign) [⌘+Shift+T]
Move one character up, down, left, or right	Arrow Keys
Move to the beginning of the line	Home
Paste a name into a formula	F3 [None]
Paste a selection	Ctrl+V [⌘+V]
Repeat the last action	F4 Or Ctrl+Y [⌘+Y]
Selects the entire worksheet	Ctrl+A
Start a formula	= (Equal Sign)
Start a new line in the same cell	Alt+Enter [⌘+Option+Enter]
Undo	Ctrl+Z [⌘+Z]