




University of Mustansiriyah
College of Pharmacy



Ali Rasool Albakaa





Laboratory Data Sheet of Experiment-1 Fundamentals of Spectrophotometer Assay of tetracycline

Name: _____ Section: _____

1. Stock solution

- (1) Mass of TC = _____ mg
- (2) Molecular weight of TC = _____ g/mol
- (3) Moles TC = mass/molar mass = _____ mol TC
- (4) Molarity of stock solution (TC) = mol/L = _____ mol/L
- (5) ppm of stock solution = _____ mg/L

2. Preparation of standard solutions

Standard solution (dilute 1.0 mL stock to 25mL):-

ppm = $(1.0/25) \times (\text{stock solution concentration}) =$ _____ mg/L

Same for concentration of solutions No.2, No.3 and No.4 with 2, 3 & 4mL of stock solution

ppm for standard solution No.2 = _____ mg/L

ppm for standard solution No.3 = _____ mg/L

ppm for standard solution No.4 = _____ mg/L

3. The absorption spectrum of TC:-

No. of record	Wavelength	Absorbance
1	400	
2	420	
3	440	
4		
5		
6		
7		
8		
9		
10		
11		
12		
13		
14		
Last	600	

Table2.for absorption spectrum of TC


Make a graph of **A** vs λ at 20 nm wavelength intervals from 400 nm to 600 nm (as measured by the whole laboratory group), at absorbance begin larger reduce wavelength intervals to 5nm, Attach the graph with this lab report, and find λ_{max} from graph?

1. The calibration curve at $\lambda_{\text{max}} =$ _____ nm

**College of Pharmacy-University of Mustansiriyah - Fifth Year Practical
advance pharmaceutical analysis – 2019 -**



Standard	Concentration	Absorbance
1		
2		
3		
4		
unknown		



On the graph you can see that the trend line does not pass exactly through zero, this is as expected: it may be due to the statistics of the data point, which are not exactly on the straight line due to random errors in the concentration and/or the absorbance reading, or to the fact that there is a remaining solution absorbance (relative to the blank) for the standards, for these reasons, you should not choose the “set intercept” box when formatting the trend line.

Finally, note that the trend line slope is the Molar Extinction Coefficient (Molar Absorptivity)! The equation for this application should be read as $A = \epsilon l C + \text{constant}$, with the constant = intercept.

2. Calculate concentration of the unknown solution from graph of calibration curve & from the equation $A = \epsilon l C + \text{constant}$?

$$C_{\text{unknown}} = C_{\text{graph}} \times 5 \quad (\text{How?})$$

3. Calculate the weight of TC in capsule?
4. Write an equation reaction (TC + Fe^{+3})?

Name: _____ Section: _____

4



1. The standard addition curve at $\lambda_{\max} =$ nm:-

Solution	Concentration	Absorbance
#1		
#2		
#3		
#4		
#5		
#6		




Figure 2: The standard addition curve

On the graph you can see that the trend line does not pass through zero, this is as expected: it may be due to the unknown addition.

2. Calculate concentration of the unknown solution from graph of a standard addition curve?

$$C_{\text{unknown}} = C_{\text{graph}} \times 10 \quad (\text{How?})$$

3. Calculate the weight of Aspirin in tablet?

4. Write an equation reaction (sodium salicylate + Fe^{+3})?



**Laboratory Data Sheet of Experiment-3
Titration of the ascorbic acid (vitamin C) in tablets
By pH meter used first and 2nd derivatives**

Name: _____ Section: _____

1. Ascorbic acid solution

(1) Mass of Ascorbic acid tablet = _____ mg

(2) Molecular weight of Ascorbic acid = _____ g/mol

(3) Normality of stock solution (Ascorbic acid) = Eq./L = _____ Eq./L

(4) ppm of Ascorbic acid solution = _____ mg/L

(5) Moles Ascorbic acid = mass/M.Wt. = _____ mol

(6) Calculate the Ka of Ascorbic acid

(7) Calculate Wt.% of Ascorbic acid

(8) 1mL of (0.1N) NaOH = ----- mg of Ascorbic acid.

No	V. of NaOH	pH	DV= v2-v1	DpH=pH2-pH1	DpH/DV	V.A . =V1+V2 / 2
1.	0.0mL		-----	----	---	---
2.						
3.						
4.						
5.						
6.						
7.						
8.						
9.						
10.						
11.						

College of Pharmacy-University of Mustansiriyah - Fifth Year Practical
advance pharmaceutical analysis – 2019 -



Using a computer data-sheet, **plot** the titration curve obtained and determine as exactly as possible the volume of titrant corresponding to the inflection point.

For more aspiring students:- The preciseness of determination of the inflection point can be very improved if calculating the second derivative of the titration curve.

Record all the numbers obtained, as well as the calculations made.

Note your observation concerning the comparison of your result with the factory value.

No	V. A.	DpH/DV	DVA= VA2-VA1	D ² pH= (DpH/DV) ² -(DpH/DV) ₁	D ² pH/DV ²	V.A.A= (VA1+VA2 / 2)
1.			-----	-----	-----	-----
2.						
3.						
4.						
5.						
6.						
7.						
8.						
9.						
10.						
11.						



Laboratory Data Sheet of Experiment-5

Titration of the Sodium Citrate in tablets

By pH meter used first and 2nd derivatives

Name: _____ Section: _____

1. Sodium Citrate solution

- (1) Mass of Sodium Citrate tablet = _____ mg
- (2) Molecular weight of Sodium Citrate = _____ g/mol
- (3) Normality of Sodium Citrate solution = Eq./L = _____ Eq./L
- (4) ppm of Sodium Citrate solution = _____ mg/L
- (5) Moles Sodium Citrate = mass/M.Wt. = _____ mol
- (6) Calculate the K_a of Sodium Citrate.
- (7) Calculate Wt.% of Sodium Citrate.
- (8) 1mL of (0.1N) HCl = ----- mg of Sodium Citrae.

No	V. of NaOH	pH	DV= v2-v1	DpH=pH2-pH1	DpH/DV	V.A . =V1+V2 / 2
1.	0.0mL		---	----	-----	-----
2.						
3.						
4.						
5.						
6.						
7.						
8.						
9.						
10.						
11.						

**College of Pharmacy-University of Mustansiriyah - Fifth Year Practical
advance pharmaceutical analysis – 2019 -**



Using a computer data-sheet, **plot** the titration curve obtained and determine as exactly as possible the volume of titrant corresponding to the inflection point.

For more aspiring students:- The preciseness of determination of the inflection point can be very improved if calculating the second derivative of the titration curve.

Record all the numbers obtained, as well as the calculations made.

Note your observation concerning the comparison of your result with the factory value.

No	V. A.	DpH/DV	DVA= VA2-VA1	D ² pH=(DpH/DV) ² - (DpH/DV) ¹	D ² pH/DV ²	V.A.A= (VA1+VA2 / 2)
1.			-----	-----	----	-----
2.						
3.						
4.						
5.						
6.						
7.						
8.						
9.						
10.						
11.						

By pH meter used first and 2nd derivatives

(8) 1mL of (0.1N) HCl = ----- mg of Sodium Salicylate.

10

**College of Pharmacy-University of Mustansiriyah - Fifth Year Practical
advance pharmaceutical analysis – 2019 -**



Using a computer data-sheet, **plot** the titration curve obtained and determine as exactly as possible the volume of titrant corresponding to the inflection point.

For more aspiring students:- The preciseness of determination of the inflection point can be very improved if calculating the second derivative of the titration curve.

Record all the numbers obtained, as well as the calculations made.

Note your observation concerning the comparison of your result with the factory value.

No	V. A.	DpH/DV	DVA= VA2-VA1	D ² pH=(DpH/DV) ² - (DpH/DV) ¹	D ² pH/DV ²	V.A.A= (VA1+VA2 / 2)
1.			-----	-----	----	-----
2.						
3.						
4.						
5.						
6.						
7.						
8.						
9.						
10.						
11.						



Experiment- 6

Infrared Spectroscopy Experiment

Name: _____ Section: _____

There are two general regions in the infrared spectrum, namely:-

(a) **Group frequency region:-**

Having a wavenumber from -----

Here, the stretching and bending vibration bonds associated with -----
or ----- are observed frequently.

(b) **Fingerprint region:-**

Having a wavenumber from-----

Here, the vibration modes depend solely and strongly on -----.

Hooke's Law: The vibration frequency of a bond is expected to increase when the
when the reduced-----decreases and also -----
increases.

Hooke's Law is expressed as:
$$\nu = \frac{1}{2\pi} \sqrt{\frac{K}{\mu}} \quad \& \quad \mu = \frac{m_1 m_2}{m_1 + m_2}$$

 ν -----, K ----- of the bond, m_1 and m_2 = -----
-----, μ -----.

Calculate m_C = mass of the carbon atom= -----

Calculate m_H = mass of the hydrogen atom= -----

Calculate m_O = mass of the oxygen atom= -----

Calculate m_N = mass of the nitrogen atom= -----

Infrared spectroscopy measures-----of IR light absorbed by
a sample and-----of the absorptions.

The vibration **frequencies** depend on

- I. The **nature of the vibration** ----- & -----).
- II. **Bond strengths.**
- III. The **masses of the atoms involved in the vibration.**

The **intensities** depend on

- I. The **change in dipole moment** that accompanies the vibration.
- II. The **number of bonds involved.**

The heavier atoms involved the lower energy.

What is the relating between ν , λ and λ with mass of atom?