

# Powder density

---

**Lab4 Industrial pharmacy**


# Powder bulk density (method 1 USP)

It simply expresses the amount (usually weight or mass) [100 g] of untapped powder in a specified interparticulate void volume (150-250 ml).

$$\rho_0 = \frac{\text{mass}}{\text{volume}} \quad (\text{unit in g/ml} = 1000 \text{ kg/m}^3)$$



As international unit (because the measurements are made using cylinders. It may also be expressed in grams per cubic centimeter ( $\text{g/cm}^3$ )).



# Composition of powders and their relation with volume

- Powders are composed of: **particles and voids.**
- The volume occupied by a given number of particles depends on how closely they are packed.
- The packing of particles depends on:
  - A- **shape & size**
  - B- **cohesiveness**
  - C- **short-range motion**
  - D- **external forces.**

# Tapped density (method 1 USP)

Describe the bulk-density of a powder (or granules) after consolidation/compression ("tapping") the container (100 or 250 ml) of powder (1 min) a measured number of times (either 300 or 250 tap), usually from a predetermined height [14 or 3 mm] .

The method of "tapping" is best described as lifting and dropping.

# Measuring Bulk density

Is determined by pouring pre-sieved ( $\geq 1$  mm) bulk drug into a graduated cylinder (250 ml with 2 ml readable) via a large funnel and measuring the volume and weight.

# Measuring tapped density

- Tapped density measured by tapping graduated glass measuring cylinder manually or by mechanical device.

**Note:** Devices that rotate the cylinder or vessel during tapping may be preferred to minimize any possible separation of the mass during tapping down.

# Types of tapping:

1. Manual tapping
2. Mechanical tapping

## 1. **Manual tapping:**

The raising and lowering of the cylinder by hand is done either without reference to the height traversed and arbitrary acceleration in both upward and downward directions; the hand remaining in contact with the cylinder at all times (hand tapping).



**(The cylinder containing the powder is tapped by repeatedly striking its base down onto a hard surface).**

## 2. Mechanical tapping

**Tap density analyzers (tap density testers)** Use an electric motor to turn a cam under a specially constructed cylinder holder. The holder secures the cylinder to a vertical shaft which runs in a low friction bearing.

**The tapping rate is expressed in taps per minute; the rate being typically a few hundred.**

The actual rate is determined by the rotational speed of the cam under the shaft/platform.



- **Digital or electromechanical counters** are usually incorporated in the device to automatically stop the cam rotation after a predetermined (yet adjustable) number of taps.

(Carry out 10, 500, and 1250 taps on the same powder sample and read the corresponding volumes  $V_{10}$ ,  $V_{500}$ , and  $V_{1250}$  to the nearest graduated unit. If the difference between  $V_{500}$  and  $V_{1250}$  is  $\leq 2$  mL,  $V_{1250}$  is the tapped volume. If the difference between  $V_{500}$  and  $V_{1250} > 2$  mL, repeat in increments such as 1250 taps, until the difference between measurements  $\leq 2$  ml).

- The **height through which the container falls** is (drop height or stroke). It is set by the distance between the highest point on the cam and the striking surface.



# How to correct the density problem:

1. Milling
2. Slugging
3. Formulation

# Benefits of bulk density:

- Important in the consideration the **size of high dose capsule** product.
- The **homogeneity of a low dose formulation** in which there are large differences in drug and excipient densities.
- knowing the dose and formulation density to determine the appropriate size for a capsule formulation.

# Tapped density and powder flowability

- The change in tapped powder volume (**related to flow properties of powders**)
- It can be determined by compressibility which is computed from powder density using the following equation:

$$\% \text{ compressibility (carr's index)} = (\rho_t - \rho_0 / \rho_t) \times 100$$

$$\text{Hausner ratio} = \rho_t / \rho_0$$

$\rho_t$ - tapped density

$\rho_0$ - bulk density

| <b>% compressibility<br/>(Carr's index)</b> | <b>Flowability</b>        | <b>Hausner ratio</b> |
|---|---------------------------|----------------------|
| <b>≤10</b>                                  | <b>Excellent</b>          | <b>1-1.11</b>        |
| <b>11-15</b>                                | <b>Good</b>               | <b>1.12-1.18</b>     |
| <b>16-20</b>                                | <b>Fair</b>               | <b>1.19-1.25</b>     |
| <b>21-25</b>                                | <b>passable</b>           | <b>1.26-1.34</b>     |
| <b>26-31</b>                                | <b>Poor</b>               | <b>1.35-1.45</b>     |
| <b>32-37</b>                                | <b>Very poor</b>          | <b>1.46-1.59</b>     |
| <b>More than 38</b>                         | <b>Extremely<br/>poor</b> | <b>&gt; 1.6</b>      |

Please follow these links to watch the videos:

- <https://www.youtube.com/watch?v=baFHheiG3Hc>
- <https://www.youtube.com/watch?v=-pEa8U4oPbA>
- <https://www.youtube.com/watch?v=72jUIIdhFcA>
- <https://www.youtube.com/watch?v=gDvEbkuhuvc>

  
**THANKS Y' ALL!**

\*\*\*\*\* YOU ARE SO KIND! \*\*\*\*\*

 **HECK YEAH!** 

**THANK YOU**

**YOU'RE THE BEST**  
*and*  
**ONLY** 