

Hydrolysis of Acetyl Salicylic Acid Solution in Sorenson Phosphate Buffer at pH 8

LAB 4

BY:


LECT. ZAHRAA AMER

LECT NORA ZAWAR

Introduction

Aspirin is a weak acid .it is soluble at 20 C°in 300 parts of water .

It is unstable in aqueous solutions degrading to salicylic acid and acetic acid



Aspirin also degrades in solid dosage forms when exposed to moisture and therefore , should be stored in tightly closed containers and kept in a dry environment .



Aim of Experiment

To study the effect of temperature on the hydrolysis of aspirin (S.P.B), and to calculate the shelf life of aspirin .

Shelf life of any drug :it is the time required for the drug to lose 10% of its effectiveness

Procedure :

1. Dissolve 0.695 g of aspirin in 250 ml of phosphate buffer (use a volumetric flask)

2. place 200 ml of this solution in an erlenmeyer flask , then keep the flask in a water bath for (30 min) at required temperature. the temperature that will be used are(40, 60, 80 °C)

Withdraw (1 ml) sample at the end of 30 mins , then continue withdrawing (1 ml) sample at 15 min interval for 90 min.

Add to each sample , 5 ml color developing reagent and read the absorbance at 530 nm

Tabulate your data.

Notes:

1- absorbance should increase with time

2- concentration of S.A at zero time is zero because aspirin not hydrolyzed yet , while the concentration of aspirin at zero time C_0 = initial concentration of aspirin 2.78 mg /ml

$$C_0 = 0.695 / 250 \times 1000 = 2.78 \text{ mg/ml}$$


3- S .A give violet color with color developing reagent

4- rate of aspirin hydrolysis follows **first order reaction**

$$\frac{dc}{dt} \propto C$$

$-\frac{dc}{dt} \propto C$ (since conc. Of aspirin decrease with time as hydrolyzed)

$$-\frac{dc}{dt} = KC$$



$$\int_{c_0}^{c_t} \frac{dc}{c} = -k \int_0^t dt$$

$$-(\ln c_t - \ln c_0) = -kt$$



$$\ln c_t - \ln c_0 = -kt$$

$$\ln c_t = \ln c_0 - kt$$

$$\text{Since } \ln = \log \times 2.303$$

$$\log c_t \times 2.303 = \log c_0 \times 2.303 - kt$$

$$\log c_t = \log c_0 - \frac{kt}{2.303}$$





5- Sorenson phosphate buffer (pH 8)

Consist of two solutions

A-1/15 M Monopotasium phosphate KH_2PO_4

B- 1/15 M Disodium phosphate $\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$

Time	Absorbance	Conc. Of S.A	Conc. Of aspirin hydrolyzed	Conc.of aspirin remained
0		$Y = c + bx$	*180/138	2.7- C_h
30				
45				
60				

7- use the calibration curve of S .A (exp. 1) to find the conc. Of S .A

$$Y = C + bx$$

Ab= intercept + slope \times conc.

e.g. $0.65 = 0.016 + 1.83 \times \text{conc.}$ $\text{Conc.} = 0.35 \text{ mg /ml of S.A}$

Conc. Of ASP. Degrade = $0.35 \times 180 \text{M. Wt of Asp} / 138 \text{ M.Wt of S.A}$

Since each mole of aspirin give 1 mole of S.A and 1 mole of A.A.

Conc. Of remaining aspirin = Initial conc. Of ASP — conc. of ASP degrade

$$2.78 - 0.46 = 2.32 \text{ mg /ml}$$

Take the log or ln conc. Remaining and plot against time in min. at each temp.

8- from the plots find the slopes then K at each temp.

9- take $\log K$ and plot $1/T$

(draw Arrhenius plot) to find K at $25\text{ }^\circ\text{C}$.


Arrhenius equation :

$$\log K = \log A - (E_a)/(2.303 R) * 1/T$$

Where A = frequency factor , E_a = energy of activation , T = absolute temp.

R = gas constant .

$$t_{10\%} = 0.105/ K_{25^\circ\text{C}}$$





THANK YOU