

Lab.4

Hydrolysis of Acetyl Salicylic Acid Solution in Sorenson Phosphate Buffer at pH 8

Presented By:

Assist. Prof. Anas Tareq

Lecturer Zeina Dawood

Lecturer Nora Zawar

Lecturer Sura Zuhair




Introduction

Aspirin is a weak acid. It is soluble at 20 C° in 300 parts of water.

Stability:


1. It is unstable in **aqueous solutions**, degrading to salicylic acid and acetic acid



2. Aspirin degrades in **solid dosage forms** when exposed to moisture [should be stored in tightly closed containers and kept in a dry environment] .



**The hydrolysis reaction is critical
It directly impacts the drug's effectiveness and shelf life.
shelf life.**



Aim of Experiment

To study the **effect of temperature** on the hydrolysis of **aspirin**, and to calculate the **shelf life of aspirin**.



Shelf life:

It is the **time** required for the drug to **lose** 10% of its effectiveness.



Procedure :




- ▶ 1. Dissolve 0.695 g of aspirin in 250 ml of phosphate buffer (use a volumetric flask)
- ▶ 2. place 200 ml of this solution in an erlenmeyer flask , then keep the flask in a water bath for (30 min) at required temperature. the temperature that will be used are(40, 60, 80 °C)
- ▶ Withdraw (5 ml) sample at the end of 30 mins , then continue withdrawing (5 ml) sample at 15 min interval for 90 min.



Notes:

- ❖ Absorbance increases with time.
- ❖ Concentration of S.A at **zero** time is **ZERO** since ASA is not hydrolyzed yet, however, the initial concentration of ASA is 2.78 mg/ml ??

Concept



Is the initial conc. of ASA 2.78 mg/ml?

- We dissolved 0.695 **g** of Aspirin in 250 ml


Therefore :

$$\frac{0.695 \text{ gram}}{250 \text{ ml}} = 0.00278 \text{ gram per ml}$$

Since 1 gm = 1000 mg

So :

$$0.00278 * 1000 = \mathbf{2.78 \text{ mg / ml}}$$





Sorenson phosphate buffer (pH 8)

Consists of two solutions:



A- 1/15M Monopotasium phosphate (KH₂PO₄)

B- 1/15M Disodium phosphate Na₂HPO₄·2H₂O

❖ The rate of ASA hydrolysis follows a **first-order reaction**



- **Use the calibration curve equation**
- **Y= C + bX (Put the equation for others)**



We are seeking ASA, not SA

Since each MOLE of aspirin give 1 mole of salicylic acid and 1 mole of acetic acid


$$\text{Mole} = \frac{Wt}{M Wt}$$

Therefore

$$\frac{wt \text{ of SA}}{M wt} = \frac{Wt \text{ of ASA}}{M wt}$$

$$\text{So : } 0.35 \text{ (from the example) } * 180 = X * 138$$

$$X = 0.456 \approx 0.46$$




- Calculate remaining ASA from:
Initial conc. of ASA – conc. of ASA degraded
 $2.78 - 0.46 = 2.38 \text{ mg}$

Time	Absorbance	Conc. of S. A	Conc. of aspirin hydrolyzed	Conc. Of aspirin remined (ct)	Log ct
0	0	0	0		
30		$Y = c + bx$	$X * 180/138$	$2.78 - Ct$	
45					
60					
75					
90					



First order Kinetics

- ▶ $\frac{dc}{dt} \propto C$
- ▶ $-\frac{dc}{dt} \propto C$ (since conc. Of aspirin decrease with time as hydrolyzed)
- ▶ $-\frac{dc}{dt} = KC$





Data Analysis: Plotting and Calculating:

Calculating:

- ▶ $\int_{c_0}^{c_t} \frac{dc}{c} = -k \int_0^t dt$
- ▶ $-(\ln C_t - \ln C_0) = kt(t-0)$
- ▶ $\ln C_t - \ln C_0 = -kt$
- ▶ $\ln C_t = \ln C_0 - kt$ ★
- ▶ Since $\ln = \log \times 2.303$
- ▶ $\log C_t * 2.303 = \log C_0 * 2.303 - kt$
- ▶ $\log C_t = \log C_0 - kt/2.303$ ★

$$Y = C + bX$$

$$\log C_t = \log C_0 - kt/2.303$$



What to do with data ?

1-take the Log, or Ln, for each calculated conc.

2-Plot the log/ Ln of conc. against time in min for each temp.

You will have 3 plots

3-Calculate K for each temp.

4-Take log K and plot $1/T$

5-(draw Arrhenius plot) to find K at 25 °C


▶ Arrhenius equation :

$$\text{Log } K = \text{log}A - (E_a)/(2.303 R) * 1/T$$

Annotations: $\text{Log } K$ is labeled 'Y', $\text{log}A$ is labeled 'Intercept', $(E_a)/(2.303 R)$ is labeled 'Slope', and $1/T$ is labeled 'x'.

Where:

- ▶ A= frequency factor,
- ▶ E_a = energy of activation,
- ▶ T= absolute temp.(temp.+273)
- ▶ R= gas constant .
- ▶ **$t_{10\%} = 0.105/ K_{25^\circ\text{C}}$**



Report to have done according to these results

- 1- Calculate K for each temp.
- 2- Calculate the Aspirin shelf life

The results for ASA hydrolysis in pH 8 was as the follow

Time in min	Absorbance at 40 C	Absorbance at 60 C	Absorbance at 80 C
0	0	0	0
30	0.095	0.125	0.33
45	0.118	0.2	0.5
60	0.15	0.27	0.64
75	0.2	0.31	0.76
90	0.23	0.4	0.76

Knowing that the calibration curve was $y = x + 0.0897$

Note: the calibration curve was in mg/ml

