Ministry of Higher Education and
Scientific Research
Al-Mustansiriyah University
College of Engineering
Environmental Engineering
Department

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

Objective:

In this experiment, the force generated by a jet of water as it strikes a flat plate or hemispherical cup may be measured and compared with moment flow rate in the jet.

Equipments and apparatus:

- The unit essentially consists of:
 - Base Plate.
 - Inlet connection.
 - Drain connection.
 - Perspex vessel.
 - Nozzle.
 - Deflector.
 - Lever mechanism.
 - Loading weights.
 - Various deflectors can be fitted at position .
 - Plate.
 - Hemisphere.
 - Cone.

Procedures:

- 1- Place the test set-up on the device so that the drain routes the water into the channel.
- 2- Fit connecting hose between device and unit.
- 3- Open device drain.
- 4- Assemble deflector , (Plate, Hemisphere or Cone). Loosen the 3 screws on the cover and remove cover together with lever mechanism. Fit appropriate deflector. Do not forget to tighten lock nut on rod. Screw cover back onto vessel.
- 5- Use adjusting screw to set pointer to zero (zero notch [7]). When doing so, do not place any loading weights on measurement system.
- 6- Apply desired loading weight 0.2N; 0.3N; 1N; 2N; 5N or combinations thereof.
- 7- Close main device cock.
- 8- Switch on device pump.
- 9- Carefully open main cock until pointer is on zero again.
- 10- Close device drain cock.
- 11- Determine volumetric flow. This involves recording time t required to fill up the volumetric tank of the device from 20 to 30 litres.
- 12- Add loading weights and note down time t for 10 litres.
- 13- Switch off pump, open drain.

Theoretical background and Calculation:

The momentum equation based on Newton's 2nd law of motion states that the algebraic sum of external forces applied to control volume of fluid in any direction equal to the rate of change of momentum in that direction. The external forces include the component of the weight of the fluid and of the forces exerted externally upon the boundary surface of control volume. If a vertical water jet moving with velocity 'V' made to strike a target (Vane) which is free, to move in vertical direction, force will be exerted on the target by the impact of jet. Applying momentum equation in X- direction, force exerted by the jet on the vane is given by:

$$F_j = \rho Q (V_{x2} \cos \beta - V_{x1})$$

The force F on the vane is equal and opposite to this, namely

$$F_j = \rho Q (V_{x1} - V_{x2} \cos \beta)$$

For flat plate, $\beta = 90$, $\cos 90 = 0$,

$$F_j = \rho Q (V_1 - V_2 (0)),$$

$$F_j = \rho Q V_1 = \rho Q V.$$

For hemispherical plate, $\beta = 180$, $\cos 180 = -1$,

$F_j = \rho Q [V_1 - (-1) V_2]$, where $V_1 = V_2$ (If we neglect the effect of change of elevation on jet speed, and the loss of speed due to friction over the surface of the vane, then $V_1 = V_2$), so

$$F_j = 2 \rho Q V$$

is the maximum possible value of force on the hemispherical cup. This is just twice the force on the flat plate.

The mass flow rate m in the jet is found by timing the collection of a known mass of water. The velocity V_1 of the jet as it leaves the nozzle is found from the volumetric flow rate Q and the cross sectional area A of the nozzle. The velocity V_2 with which the jet strikes the vane is slightly less than V_1 because of the deceleration due to gravity. This effect may be calculated from the expression:

$$V_2^2 = V_1^2 - 2gs$$


Inserting the value $s = 0.035$ m leads to the result

$$V_2 = \sqrt{V_1^2 - 2gs}$$


$$F_{\text{mass}} = \text{mass} * g (9.81 \text{ m/ sec}^2)$$

Where $Q =$ Discharge from the nozzle (Calculated by volumetric method)

$V =$ Velocity of jet = (Q/A)

Item	Value
Nozzle diameter:	8mm
Distance between nozzle & target plate:	
Diameter of target plate:	36mm
Target plates:	- 180° hemispherical target - 120° target (cone) - flat target - 30° target

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