# Experiment NO. 3

### Determination of Nickel as Dimethyl glyoxime complex.

#### **Introduction**

It is considered one of the most older Organic precipitation and The optimum of them.

It is a weak acid  $(CH3C = NOH)_2$  and dissolves slightly in water and gives one ion of hydrogen when ionized:

This Detector is dissolved in alcohol so alcohol solution is used Specially to precipitate nickel quantitatively.

Ni<sup>+2</sup> is precipitated by unites with two molecules of DMG and releasing two of hydrogen ions that is equivalent with an increase of ammonium hydroxide according to the following reaction:

$$2 \begin{array}{c} \text{CH}_{3}\text{-C=N-OH} \\ 2 \begin{array}{c} | \\ \text{CH}_{3}\text{-C=N-OH} \end{array} + \text{Ni}^{+2} + \text{NH}_{4}\text{OH} \longrightarrow (\text{C}_{4}\text{H}_{7}\text{O}_{2}\text{N}_{2})_{2}\text{Ni} + 2\text{NH}_{4}^{+} + 2\text{H}_{2}\text{O} \end{array}$$

The precipitate is red color and its solubility in water very low. we can dry it in  $(110-120 \, ^{\circ} \, \text{C})$  but is dissolves in dilute mineral acids .

Hydrogen released at the Reagent Union with nickel due to increased solubility of precipitate .

## Materials

- 1- Sample containing nickel.
- 2 Dilute ammonia solution.
- 3 Dilute hydrochloric acid.

#### **Procedure**

- **1-** Carefully weigh (0.1 g) of Nickel in the beaker of 400 ml with a glass rod and Dissolves in the least amount of distilled water and stirring the solution To the dissolve.
- **2-** Add (2 mL) of dilute hydrochloric acid (1: 1) and add distilled water to 75 ml.
- **3-** Heat the solution on the heater and add (25 )ml of the precipitate agent (DMG) hot solution. Add the dilute ammonia solution quickly as a dropes with a continuous stirring until is fully precipitation.
- **4-** Put the beaker with its content on the heater for a quarter of an hour until a red precipitate appears.
- **5-** Prepare filtration device, and weigh the empty filter paper, then filtered the solution and wash it with cold water.
- **6-** Dry the paper with the precipitate and then weigh it and calculate the percentage of Nickel.

## Calculation

G.F = 
$$\frac{(A.Wt.(Ni))}{(M.Wt.Ni(DMG)2)} = 0.20314$$

$$Ni\% = \frac{(wt.Ni(DMG)2 \times G.F)}{(wt.of sample)} \times 100$$

## **Questions of discussion**

- 1- What are DMG properties? And why it is used as a precipitating factor?
- 2- What are Ni (DMG) <sub>2</sub> specifications?
- 3- Add dilute hydrochloric acid?
- 4- Why was a precipitate process of Ni (DMG) 2 in the base medium?
- 5- What is the chemical composition of the Ni (DMG) 2 precipitate?