

Experiment No.(1):- Preparation and standardization of 0.1 M(HCl) hydrochloric acid solution

Theory:- Hydrochloric acid is produced in solutions up to 38% HCl (concentrated grade). Higher concentrations up to just over 40% are chemically possible, but the [evaporation](#) rate is then so high that storage and handling need extra precautions, such as pressure and low temperature. Laboratory grade hydrochloric acid is not sufficiently pure to be used as a primary standard, because it evaporates easily. In this experiment, a standard solution of sodium carbonate is used to determine the exact concentration of a hydrochloric acid solution. The neutralization reaction that occurs is as follows:



Methyl orange indicator solution is used. At the end-point – when neutralization just occurs – the indicator changes color from yellow to peach-pink.

Procedure:-

1. [Preparing \(50 ml\) 0.1 M HCl Solution](#): 38 % HCL shows density 1.19 g/mL and we can find M by next :-

$$M = \frac{\text{sp.gr} * \% * 1000}{\text{M.wt}}$$

Calculate the volume of HCl (conc.):- We must dilute it to preparing 0.1 M HCl in 50 ml from next:

$$(M * V)_{\text{conc.}} = (M * V)$$

dilute

$$M * V_{\text{ml}} = 0.1 * 50 \text{ ml}$$

Transfer V ml by cylinder to clean and dry beaker containing 30 ml D.W, transfer the solution to volumetric flask capacity 50 ml, and complete the volume to the mark by D.W.

2. [Preparing \(50 ml\) 0.1 M Na₂CO₃ Solution](#):-calculate amount from sodium carbonate for prepare 0.1 M in 50 ml -

$$M = \frac{\text{Wt. (gm)}}{\text{M.Wt.}} * \frac{1000}{V(\text{ml})}$$

$$0.1 = \frac{\text{Wt. (gm)}}{106} * \frac{1000}{50}$$

$$\text{Wt.} = 0.53 \text{ gm}$$

Weigh 0.53 gm. from Na_2CO_3 in clean and dry beaker and dilute in 30 ml D.W, transfer solution to volumetric flask capacity 50 ml and complete the volume to the mark by D.W.

3. Transfer known volume of 5 ml the sodium carbonate solution, with a pipette, to a conical flask then add one or two drops of methyl orange to this solution.

4. Add the acid unknown solution from the burette gradually with continuous swirling of the solution in the conical flask and near the end point, the acid is added drop by drop. Continue the addition of the acid until the color of the solution passes from yellow to faint red.

5. Repeat the experiment three times and tabulate your results then take the mean of the three readings.

Calculations: Calculate the molarity of HCl:-

$$m \text{ mol HCl} = m \text{ mol Na}_2\text{CO}_3$$


$$(M * V) \text{ HCl} = (M * V) \text{ Na}_2\text{CO}_3 * \frac{1}{2}$$

$$(M * V_{\text{burette}}) = (0.1 * 5) * \frac{1}{2}$$

Discussion:-

1. What the difference between primary and secondary standard substances?

2. Calculate the volume of conc. HCl required for preparing 250 ml 0.1 M?



3. Calculate the weight of Na_2CO_3 required for preparing 100 ml 0.1 M?
4. Why is sodium carbonate primary solution?
5. Why standard solution should be colorless?
6. Why is HCl not primary solution?
7. What is the titration?