**Al-Mustansiriyah University**  
**College of Science/Department of Chemistry**

**Analytical Chemistry Lab.**  
**Second Year**  
**Four Lecture Exp.2 /II Semester**

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

## Experiment (2)

- Determination of the number of total equivalents of salts
- Determination the ratio of the NaCl and KCl salts using cationic exchanger

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
### Materials

- 1- HCl (3M)
- 2- AgNO<sub>3</sub> (0.1 M)
- 3- Mixture Saline (NaCl + KCl)
- 4- Methyl orange (M.O)
- 5- Standard base NaOH (0.1N)

## Procedure

1. **Activate** the **cationic column** by hydrochloric acid (**3M HCl**), as we know

2. Carefully **weigh (0.2 g)** from the sample contain (**NaCl & KCl**)

3. **Dissolve** the salt with a little amount of the **distilled water** in the beaker, then transfer it to the **cationic** exchanger

4. **Collect** the solution from the **bottom** of the column (eluent solution) in the **conical flask 1** and then wash the beaker with a little distilled water and transfer it to the ion exchange column.

5. Collect the first amount in **conical flask1** of washing water from the column and check it by adding of the **M.O indicator**. when a **pink color appear**, transfer to conical **flask2** Add another amount of the distilled water to the column and we take **eluent solution** it in the conical flask1, then also check it by use the M.O indicator and **repeat** this process until the **pink color disappearance** and it **changes yellow color**
6. **Titrate** the solution that **collected pink color** in the conical flask **2** with the sodium hydroxide **NaOH (0.1N)** .
7. We **find** the standard base **volume** by a **burette** Let it be (**20 ml**).

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1. No. of meq. of cation (**NaCl + KCl**) = No. of meq. NaOH

No. of meq. of NaCl + No.of meq. of KCl = No. of meq. of NaOH

$$N_{(\text{NaCl})} \times V_{(\text{NaCl})} + N_{(\text{KCl})} \times V_{(\text{KCl})} = N_{(\text{NaOH})} \times V_{(\text{NaOH})\text{from burette}}$$

$$\frac{\text{wt}_{\text{NaCl}}}{\text{eq. wt}_{\text{NaCl}}} + \frac{\text{wt}_{\text{KCl}}}{\text{eq. wt}_{\text{KCl}}} = \frac{N_{(\text{NaOH})} \times V_{(\text{NaOH})}}{1000}$$

$$\text{Wt of NaCl} = (X) \text{ g}$$

$$\text{Wt of KCl} = (0.2 - X)\text{g}$$

$$\frac{x}{58.5} + \frac{0.2 - x}{74.5} = \frac{N_{(\text{NaOH})} \times V_{(\text{NaOH})}}{1000}$$

$$2. \text{NaCl \%} = \frac{x}{0.2} \times 100$$

$$3. \text{KCl \%} = \frac{0.2 - x}{0.2} \times 100$$

## Discussion questions

- 1- How do we know that the column has become **neutral**?
- 2- What will happen to the column after we transfer the solution sample to it and what the type of the groups that the column will be carry?
- 3- What is the **eluent solution** that collected from the bottom of the **column cationic** if the component solution was passed through the resin

(CaCl<sub>2</sub>, H<sub>2</sub>SO<sub>4</sub>, H<sub>2</sub>O , KOH, F)

**4-** (0.5)g of mixture saline solution (30%silver chloride and potassium nitrate ) was passed through a cationic exchanger and the collected solution from the column was calibrated against the (0.1N) sodium hydroxide base, the equivalent size of the eluent solution (10ml), **calculate**

1. The weight of AgCl and weight of KNO<sub>3</sub> in the sample?
2. The percentage of Cl in AgCl?
3. ppm for potassium nitrate in the solution (5ml)?
4. Write the equation for this question?

The atomic weight for (K = 39 , N =14 , O = 16 , Cl = 35.5 , Ag = 108 , Na = 23)