Lab : 9

Soil Microbiology

Nitrogen Cycle

Nitrogen fixation :

This process in which nitrogen (N_2) in the atmosphere is convert into ammonia (NH_3) fixed form ,by number of soil microorganisms are capable to production an enzyme called nitrogenase which needs to metal molybdenum and iron to complete the process. This process is essential for all forms of life because nitrogen is required to biosynthesize basic building blocks of plants, animals and other life forms, e.g., nucleotides for DNA and RNA and amino acids for proteins.

Two types of bacteria are responsible for fixing nitrogen :

Family 1 : Azotobacteriaceae ex. Azotobacter.

Family 2 : Rhizobiaceae ex. *Rhizobium*.

Some types of bacteria can fixing nitrogen but in less efficient such as *Clostridium*, *Klebsiella*, Cyanobacteria.

Nitrogen fixation can possible to get in two ways:

1- Symbiotic Nitrogen fixation . ex. Rhizobium.

2- Non-Symbiotic Nitrogen fixation (free living). ex. . Azotobacter sp. , Azospirillum

Symbiotic Nitrogen fixation :--- ex. Rhizobium

Nitrogen fixation through nodule formation in Leguminous plants

Leguminous plants bearing nodules on its root as a result of injury root hairs with bacteria *Rhizobium*, In large numbers inside this root nodules. Bacteria be parasite on plant where they get the necessary energy for the growth and fixed nitrogen that useful for plant. The symbiotic Bacteria Rhizobia (from the Greek words Riza = Root and Bios = Life)

Procedure:

- 1- Take plant called sweet clover (Melilotus) from soil.
- 2- Cut roots & wash it well with tap water to remove soil particles .

Lab : 9

Soil Microbiology

- 3- Cut large pink or red nodules true nodules (neer the root), False nodules Located far from root.
- 4- Put true nodules on clean slide and smash it with another slide , leave it to dry and stain with gram pigment .
- 5- Examine under oil immersion, *Rhizobium* appear look like latin letters Y, X, Z, T.

Non-Symbiotic Nitrogen fixation (free living):

The fixation of free nitrogen in the soil can gets by most microorganisms living freely or outside the plant cell that called non-symbiotic N_2 fixation. It is performed by the aerobic and anaerobic bacteria and blue green algae.

The process involved microorganism in soil free-living *Azotobacter chroococcum* bacteria , this bacteria G-ve, motile , often pleomorphic ranging from rode to cocci shapes , strict aerobic, the colonies appear like sticky mucus because they around by cyst , when bacteria grown on agar media form large convex , mucoid colonies , with brownish color, growth in PH (7.2-8) .To isolation this bacteria use (Nitrogen –free glucose medium) ,this contains sources carbon, molybdenum .

This bacteria produce two type of pigments :Use to distinguish between 6 types of *Azotobacter*

1-Water- soluble pigments. (It's the important because it spread in culture medium).

2-Water-in soluble pigments.

Procedure:

1-Transfer 1gm of soil to Nitrogen –free glucose broth medium that packaged in glass bottles.

2- Close the bottles not tight to provide enough aeration and get N_2 from air .

3- Incubate at 30c for 7days in horizontal position to increase the surface area of culture medium that exposed to the air .

4- after incubation :transfer 0.1 ml from Nitrogen – free glucose broth to Nitrogen – free glucose agar and streaking to show colony shape and color water- soluble pigment.