

The nitration mechanism of benzene :



Chemical and Apparatus

Benzene, conc.HNO₃, conc.H₂SO₄, cold distilled water, ethanol, ice path, Round bottom, condenser, funnel, electric heater, beaker, balance, filter paper, oven,

Safety Note

Caution: Avoid contact with the acids used in this experiment and the reaction product. Prevent contact with the skin, eyes, and clothing; work in the hood. An acid spill is neutralized using solid sodium carbonate or bicarbonate. The reaction is highly exothermic.

Concentrated nitric acid and concentrated sulfuric acid are both <u>strong oxidizers</u>, and <u>strongly corrosive</u>--wear gloves while handling them, and avoid breathing their vapors.

Experimental No. (10)

Procedure of Experimental

- 1- In hood, prepare a mixture of (4 mL) conc. HNO₃ and (4 mL) conc. H₂SO₄ in a beaker. Cool beaker to room temperature (30^oC) by means of a ice bath, , remove the beaker from the ice bath, wipe off the outside of the conical flask, and clamp it to hotplate .
- 2- To the conical flask, add (0.5 mL) benzene <u>DROPWISE</u> over a period of about 5-10 min; gently swirling the tube to mix the contents after each addition. Keep the reaction mixture between 50 55 °C. <u>DO NOT ALLOW THE REACTION MIXTURE</u> <u>TO EXCEED 60 °C</u>.
- 3- After the addition is completed and the exothermic reaction has subsided, put the mixture in the round bottom with condenser and heat the round for (45 min) in a hot water bath, maintaining the temperature in the reaction below 90 °C during this period.



Figure show the reflux of reaction mixture

- 4- Cool the conical flask in an ice bath to room temperature.
- 5- Pour the reaction mixture into 100 mL of distilled water which is in a 150 mL beaker.

- 6- Isolate the crude product by filtration.
- 7- Wash the filter cake thoroughly with cold (0-10°C) distilled water and dry the filter in the oven.
- 8- The crude product may be purified by recrystallization by hot ethanol, then determine the weight, melting point, and percentage yield.

Calculation



Theoritical weight =	Wt of benzene \times M.Wt of nitrobenzene
of nitobenzene	M.Wt of benzene

Experimental weight = weight of product with filter paper – weight of filter paper

of nitrobenzene

% of nitrobenzene = **Experimental weight X** 100 **Theoritical weight**

Questions for discussion

- 1. Which is nitrated faster? toluene or nitrobenzene, toluene or phenol? Explain.
- 2. List three combinations of reagents used for nitration of aromatic compounds?
- 3. Why is concentrated sulfuric acid employed in this reaction? What is the electrophile that is produced by the reaction of sulfuric acid and nitric acid/
- 4. Why is it important to maintain the reaction temperature low and the addition of nitric acid-sulfuric acid mixture carried out slowly?
- 5. Explain why concentrated H₂SO₄, not concentrated HCl, was used in this experiment?