# Experimental No. (14)Derivatives of carbonyl compoundsDerivatives of carbonyl compounds

#### **Purpose of experimental**

Preparation of 2,4-dinitrophenylhydrazone



Aldehyde : R = alkyl or aryl, R = hydrogenKetone : R = alkyl or aryl, R = alkyl or aryl

#### Theory part of experimental

carbonyl compounds (aldehydes and ketones) react with ammonia derivatives such as 2,4-Dinitrophenyl hydrazine, semicarbazide and hydroxylamine, etc. to form orange-red coloured stable products which have well defined melting points. These compounds are useful in identification and characterization of carbonyl compounds. The reaction catalyzed by protonic acids such as HCl, H<sub>2</sub>SO<sub>4</sub>, CH<sub>3</sub>COOH etc. involves nucleophilic addition of the ammonia derivatives to the carbonyl carbon followed by elimination of water molecule. pH of the reaction should be maintained around 5 to 6.

#### **Chemical and Apparatus**

Benzaldehyde, 2,4-dinitrophenylhydrazine, ethanol, sulphuric acid, beaker, heater, filter paper, funnel,

## **Procedure of Experimental**

## Part A: Preparation the 2,4-dinitrophenylhydrazine test reagent.

- Weigh out 0.5 g of 2,4-dinitrophenylhydrazine and dissolve in 10 mL of concentrated sulphuric acid in a beaker.
- 2. Add this solution carefully to a solution of 15 mL of ethanol and 5 mL of deionised water. This gives the yellow-brown stock solution. The solution can be used for about 10 days but it gradually deteriorates as a precipitate forms

# Part B: Preparation of Benzaldehyde 2,4-Dinitrophenyl hydrazone

- 1- In a clean and dry test tube or beaker take 2 ml of benzaldehyde.Add 10 ml of 2,4-dinitrophenylhydrazine reagent
- 2- Warm the contents of tube or beaker on water bath for 5 minutes and allow the tube or beaker to stand at room temperature for 5 minutes.
- 3- Cool the contents of tube or beaker in ice water bath when an orange-coloured precipitate of Benzaldehyde 2,4-Dinitrophenyl hydrazone will separate out.
- 4- Filter the product and dry it.
- 5- Record the yield of the product and determine the melting point of the product, (lit. 237 °C).



## **Calculation** $NH-NH_2$ $NO_2$ сно $NO_2$ +CH=N-NH -NO<sub>2</sub> (II) $\dot{N}O_2$ **(I)** M.Wt =M.Wt =Wt. = Wt. = weight = density x volume Wt. of (I) X M.Wt of (II) Theoritical weight of (II) M.Wt of (II) Experimental weight of (II)

% of benzaldehyde 2,4-dinitrophenylhydrazone = Experimental weight of (II) Theoritical weight of (II) X 100

## **Questions for discussion**

- 1. Can you write equations for any reactions occurring?
- 2. What is the purpose of making derivatives of unknowns?
- 3. Why we added the sulphuric acid in preparation of Benzaldehyde 2,4-Dinitrophenyl hydrazone
- 4. Which method use to purification of yield (Benzaldehyde 2,4-Dinitrophenyl hydrazone)