EXPERIMENT NO.5

Conductometric titration (acid – base)

The purpose of experiment:

Part A: Calibration of NaOH using strong acid HCl

Part B: Determination the concentration of weak acid CH₃COOH


Part A:

Apparatus: Conductivity meter, magnetic stirrer, stirrer bar, burette, beaker.

Reagents: NaOH solution, standard solution of HCl and D.W.

Procedure

- 1- Calibrate the conductometer with standard solution
- 2- Put () ml, () N HCl was peptide and transferred to 250ml beaker, diluted with 150ml D.W., add the stirrer bar and dipped the conductivity cell, we make sure that when the stirrer is turned on, the stirrer bar will not hit the cell.
- 3- The solution was stirred for a few seconds, then the initial conductance was read.
- 4- Approximately 0.2ml of NaOH solution was added at each time from the burette to the solution, stirred for a few seconds, then the conductance was read.
- 5- The titration was continued until suitable volume (5-6)ml of NaOH was added and the conductance was read after each addition.
- 6- Tabulate all the results obtained in this titration.
- 7- Plot the titration curve between conductivity versus volume of NaOH add as atitrant.
- 8- Show and determine the equivalent point and volume of NaOH at the equivalent point
- 9- Calculate the concentration of NaOH using this equation:
$$(N \times V)_{\text{HCl}} = (N \times V)_{\text{NaOH}}$$
- 10- Remove the conductivity cell and rinsed with distilled water.

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
Reagents: CH_3COOH solution (unknown), NaOH solution and D.W.

Procedure:

- 1- A () ml of unknown acetic acid was pipetted and transfer to 250ml beaker and diluted with 150ml D.W. and put the magnetic bar and dipped the electrode, make sure that when the stirrer is turned on, the stirrer bar will not hit the conductivity cell, the solution was stirred for a few seconds, then the initial conductance was read.
- 2- Approximately 0.2ml of NaOH 0.1N solution was added at each time from a burette to the solution, stirred for a few seconds then the conductance was read
- 3- The titration was continued until suitable volume of NaOH was added and the conductance was read after each addition.
- 4- Tabulate all the results obtained in this titration
- 5- plot the titration curve between conductivity versus volume of NaOH added as titrant
- 6- Show and determine the equivalent point and volume of NaOH at the equivalent point
- 7- Calculate the concentration of CH_3COOH using this equation:

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