

Carbohydrates

- **Carbohydrates:** Are the most abundant and diverse class of organic compounds occurring in nature.
- They play **a key role in the evolution of life** due to creating a direct link between the sun and chemical energy.

The aim of the experiments:

- Understanding simple tests for **the identification of carbohydrates** in a given sample.

Theory:

- **Carbo**=Carbon, **Hydrate**=Water (Hydrogen+Oxygen)
- **General formula:** $(CH_2O)_n$
- They have important **structural and metabolic** roles in both **animals and plants**.
- Commonly used for **food and energy storage**.

Classification of carbohydrates

- **Monosaccharides:** Cannot be hydrolyzed into a more simple form.
- **Disaccharides:** Products of **two** monosaccharide units.
- **Oligosaccharides:** Products of **three to ten** monosaccharides.
- **Polysaccharides:** Products of **more than ten** monosaccharide units.

Examples of carbohydrates

- Monosaccharides:
- Can be further classified into:


No. of carbon atoms Trioses

Tetroses – Pentoses – Hexoses – Heptoses

Aldehyde or Ketone group

Aldoses – Ketoses

Monosaccharides structure:




Monosaccharides examples:

	Aldoses	Ketoses
Trioses (C ₃ H ₆ O ₃)	Glyceraldehyde	Dihydroxyacetone
Tetroses (C ₄ H ₈ O ₄)	Erythrose	Erythrulose
Pentoses (C ₅ H ₁₀ O ₅)	Ribose	Ribulose
Hexoses (C ₆ H ₁₂ O ₆)	Glucose	Fructose

Other examples of carbohydrates include:

- **Disaccharides:**
- Lactose, Maltose, and Sucrose.



- **Polysaccharides:**
 - Starch and Dextrin.

Can be linear or branched polymers.

Also, sometimes they are classified as **Hexosans** or **Pentosans** depending on their yield when hydrolyzed.

Qualitative and Quantitative tests for carbohydrates

To identify the presence of carbohydrates in a sample:

- 1 - Molisch's test
- 2 - Fehling's test
- 3 - Benedict's test
- 4 - Tollen's test
- 5 - Iodine test


Carbohydrate identification tests

1 - Molisch's test:

General test for carbohydrates.

Concentrated **sulfuric acid** converts the carbohydrate into **furfural** or any of its' derivatives, which in turn **reacts with α -naphthol** to form a **purple colored product**.

Purple or violet ring is a positive result for the presence of carbohydrates in the sample.

Molisch's chemical reaction:**Molisch's Procedure:**

1. Take 2ml of the given sample solution in a clean test tube.
2. Add 2-3 drops of Molisch reagent slowly.
3. Add concentrated sulfuric acid along the sides of the test tube.
4. Notice the **separation** of the sample (Acid: bottom, sample: top).
5. If a **violet ring** is observed in the middle, then the **presence of carbohydrates** is confirmed.


2. Fehling's test:

This test identifies **reducing sugars**.

Needs **boiling water bath**.

The copper ions in Fehling's solution (+3) state is reduced to (+2) oxidation state, and **red cuprus oxide** is precipitated in the **alkaline medium**.

Red precipitate is a positive result for the presence of reducing sugars and carbohydrates.

Fehling's chemical reaction:**Fehling's Procedure:**

1. Take 2ml of the given sample solution in a clean test tube.
2. Add 1ml of Fehling A solution slowly.
3. Add 1ml of Fehling B solution slowly.
4. Keep the solution in a boiling water bath for ≈ 10 minutes.
5. If a **red precipitate** is observed, then the **presence of carbohydrates and reducing sugars** is confirmed.


3. Benedict's test:

This test identifies **reducing sugars**.

In **Alkaline medium**, **sodium carbonate** converts glucose to **enediol**, which in turn reduces **cupric to cuprous** forming **cuprous hydroxide**.

Upon boiling, **red precipitate of cuprous oxide** is formed.

Red precipitate is a positive result for the presence of reducing sugars and carbohydrates.

Benedict's chemical reaction:

Benedict's Procedure:

1. Take 3ml of the given sample solution in a clean test tube.
2. Add 3ml of Benedict's reagent.
3. Boil the solution for ≈ 2 minutes.
4. Cool the solution and observe the test tube.
5. If a red/green/yellow precipitate is observed, then the presence of carbohydrates and reducing sugars is confirmed.


4. Tollen's test:

This test identifies **reducing sugars**.

Reacts with carbohydrates to form a **silver mirror** along the inner walls of the tube.

Silver ions are reduced to **metallic silver**.

Silver mirror is a **positive result** for the presence of reducing sugars and carbohydrates.

Tollen's chemical reaction:**Tollen's Procedure:**

1. Take 3ml of the given sample solution in a clean test tube.
2. Add 3ml of Tollen's reagent.
3. Keep the test tube in boiling water for ≈ 10 minutes.

