



Practical Microbiology



3<sup>rd</sup> grade

Second course\_Lab 3

# ISOLATION AND IDENTIFICATION OF STAPHYLOCOCCI


م.م. زهراء عباس صحن

- These bacteria are G+ve.
- cells aggregation like bunch of grapes.
- non-motile.
- non-spore-forming.
- facultative anaerobic.
- colonies are round, convex, mucous.
- Staphylococci are found on the skin mouth and upper respiratory tract. Such as pathogenic *Staphylococcus aureus*.




- It may be isolated from Skin or mucous membrane of oral cavity.
- Common skin infections caused by *Staph. aureus* include: pimples, furuncle and impetigo.
- Systemic infections causes various diseases (pneumonia, endocarditis, tonsillitis, bacteremia, osteomyelitis, and food poisoning)


## Staph Infection: Skin-Related Symptoms



Boils



Cellulitis



Infection in hair follicles



Impetigo



Paronychia

- (deep tissue) infections that result from *Staph. aureus* also produce infections of the sinuses (sinusitis) and middle ear (otitis media).
- May also cause oral infection, angular cheilitis together with Candida at the angles of the mouth. In gingival carbuncles forming abscesses.




- *Staphylococcus epidermidis* is a non-pathogenic and member normal microflora of the skin.
- It may cause serious infections, if it has unusual opportunity to enter surface barriers for example cardiac surgery patients or those with intra venous catheters.

# Identification


- 1. Under microscope**
- 2. Colonial**
- 3. Catalase tests**
- 4. Coagulase test**
- 5. Protein A latex agglutination test**
- 6. DNase test:**

# Identification

## 1. Under microscope




## 2. Colonial



## Procedure:

- Swab from the throat around the tonsillar area or rotate a moist with (saline) swab around entire perimeter of both nares.
- Roll the swab near the edge of a blood agar plate, and mannitol salt agar.
- Streak with sterile loop, label the plates.
- Repeat the same steps for Gingivae and tongue label the plate.
- Incubate the plates at 35°C for 24-48 hr.




Observe the growth on **Blood agar**:

*Staph. aureus*: colonies are large , smooth  **$\beta$ -hemolytic** (clear zone around the colonies).




*Staph. aureus* on manitol salt agar:

Colonies and culture medium turn to yellow due to fermentation of mannitol.



### 3. Catalase tests


Some bacteria contain flavoproteins that reduce O<sub>2</sub> resulting in the production of hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) or superoxide (O<sub>2</sub><sup>-</sup>). These are **toxic** for obligate aerobes and facultative anaerobes. Many bacteria produce enzymes to protect them self against superoxide (O<sub>2</sub><sup>-</sup>), these enzymes are catalase, peroxidase or superoxide dismutase, which catalyze the destruction of H<sub>2</sub>O<sub>2</sub> or O<sub>2</sub> as follows:



Catalase production can be noted by mixing H<sub>2</sub>O<sub>2</sub> with the tested bacteria. Bubbles of O<sub>2</sub> represent a positive catalase test, and the absence of bubbles represents a negative catalase test. Most strict anaerobes lack both enzymes and cannot tolerate O<sub>2</sub>

## Procedure

- Label two nutrient broths or nutrient agar slants.
- Using aseptic technique, inoculate one agar slant or broth with *staphylococcus* and the other tube with *streptococcus*.
- Incubate all the cultures at 37°C for 18 to 24 hours.
- Add few drops of 3% H<sub>2</sub>O<sub>2</sub> over the growth of the slants or broths and observe the appearance of gas bubbles.
- Record the results.





## 4. Coagulase test

- This test is used to distinguish between **pathogenic** and **non-pathogenic** *staphylococcus*
- Coagulase is enzyme produced by pathogenic *Staphylococci*, that clot blood plasma . Citrate is usually added to act as anticoagulant and prevent false positive results.

### Procedure:

- Add 0.5 ml of citrated rabbit plasma to two small test tubes; label the tubes with name of bacteria.
- Add 0.5 ml broth culture of *S. aureus* to one tube and 0.5 ml broth culture of *S. epidermidis* to the other tube.
- Incubate the broth cultures at 37°C for 1 to 4 hours in water bath.
- Examine the broth cultures for the presence (**coagulase positive**) or absence (**coagulase negative**) of clouding and clots.



## 5. Protein A latex agglutination test


Protein A, synthesized by almost all strains of *Staph. aureus*, has a special affinity to fragment of immunoglobulin G ( IgG ). Hence when latex particles coated with IgG are mixed with a suspension of organism visible agglutination of latex particles occurs.

## 6. DNase test:

- Most coagulase positive *staphylococcus* produces this enzyme.
- Most pathogenic strains of *staphylococcus* produce a nuclease enzyme called DNase degrades host DNA and increases the pathogenicity of *staphylococcus*.

## Procedure:

- Divide a DNase agar plate in half with a wax pencil and label each part with the name of bacteria .
- On one part of the agar plate, put heavily spot-inoculate of *S.aureus* and the other part inoculated by the same way with *S.epidermidis*.
- Incubate at 37°C for 18 to 24 hours.
- Flood the DNase test agar with 1 N HCl. A zone of clearing around the colony indicates a positive DNase test.
- Record the result.



# Comparison between Staph. aureus & Staph. epidermidis

## *Staph. aureus*

pathogenic

Coagulase **+ve**

On mannitol salt ager:

mannitol fermentation **+ve**

**Tolerate** High Salt concentration

DNase **+ve**

## *Staph. epidermidis*

Non - pathogenic

Coagulase **-ve**

On mannitol salt ager:

mannitol fermentation **-ve**

**Will not tolerate** High Salt  
concentration

DNase **-ve**



THANKS FOR YOUR CONSIDERATION

