


A ten-fold dilution

A ten-fold dilution, also known as a 1:10 dilution, reduces the concentration of a solution or substance by a factor of ten. This means the new concentration is one-tenth of the original concentration. It's a common technique in scientific experiments, especially when dealing with very high concentrations.

- **Ratio:** A ten-fold dilution involves a 1:10 ratio, meaning 1 part of the original solution is mixed with 9 parts of diluent (water or buffer) to create a total volume that is ten times larger.
- **Concentration Reduction:** Each ten-fold dilution reduces the concentration by a factor of ten.
- **Example:** If you start with a 1 M solution and perform a ten-fold dilution, the resulting solution will have a concentration of 0.1 M (1/10 of the original).


Serial Dilutions: When multiple ten-fold dilutions are performed in a series, it's called a **ten-fold serial dilution**. This is often used to achieve a much lower concentration from a high starting concentration in a few steps.


Serial Dilution



Science Experiment ● ● ●

SERIAL DILUTION





Understanding Basics of Dilution

Let's say you have a glass of sugar solution at 6.0 M concentration and another glass with water. If you took 1 mL of sugar solution and dispensed it in a new tube, then added 9 mL of water and mixed. You just made a dilution.

To calculate the dilution factor (df), you need to divide the solution's total volume by the sample volume (i.e., sugar solution in this case).

Total volume = (Volume of sample + Volume of diluent),
i.e., 1 mL sugar solution + 9 mL water = 10 mL

$$\text{Dilution factor} = \frac{\text{Total volume}}{\text{Sample volume}} = \frac{10 \text{ mL}}{1 \text{ mL}} = 10$$

Dilution Factor = Total volume/Sample volume = 10/1= 10

This dilution can be expressed by various terms like

- **The dilution factor is 10 (dilution factor is the reciprocal of the dilution)**
- **It was a 10-fold dilution**
- **It was diluted by 1/10.**

Calculating the concentration

To calculate the concentration of this diluted sample, multiply by the inverse of the dilution factor.

Final concentration = (Initial concentration X 1/ dilution factor), i.e., 0.6 M

$$\begin{aligned}\text{Final concentration} &= \text{Initial concentration} \times \frac{1}{\text{Dilution factor}} \\ &= 6\text{M} \times \frac{1}{10} \\ &= 0.6\text{M}\end{aligned}$$

Sometimes, you may wish to calculate backward. Let's say you have a sample that has been diluted 1/10 with a concentration of 0.6 M and wants to calculate its undiluted concentration.

Undiluted concentration = (Diluted concentration X dilution factor)

$$= 0.6 \text{ M} * 10 = 6\text{M}.$$

Multiple or Serial Dilution

Let's say you again diluted the sugar solution by taking 1 ml of diluted sugar solution in another 9 mL of water. Your final dilution ends up being $1/10 \times 1/10 = 1/100$ dilution. In other words, your sample has been diluted 100-fold (d.f. =100).

If you want to measure the concentration of the final solution after two successive 10-fold dilutions, follow the same formula: $(6 \text{ M} \times 1/10 \times 1/10) = 6/100 = 0.06 \text{ M}$

So, If you need 100-fold (10^{-2}) dilution, you can either add **0.1 mL** sample with **9.9 mL** of diluent or make two successive 10-fold dilutions.

Preparation of diluent

- 1. Prepare six test tubes that can store 20 mL or more in a rack and label them T1-T6.**
- 2. Each tube is consistent with the dilution factor it corresponds to (i.e., T3= 1×10^3 OR 0.001th of stock concentration).**
- 3. Pipet 9 mL of sterile water, saline, or broth into each of the 6 test tubes.**
- 4. Sterilize the diluent banks by autoclave. Use aluminum foil to cover each of the six test tubes and then transfer them to an autoclave-compatible test tube rack. Sterilize for a minimum of 15 minutes at 121°C, 15 PSI.**
- 5. Remove blanks using heat-resistant gloves and allow to cool. When tubes reach room temperature, cover and store at 4°C until needed.**

Serial Dilution Method

- 1. Obtain the sample flask from the incubator and shake it vigorously.**
- 2. Pipet 1 mL of “solution” into the test tube labeled T1. Vortex T1.**
- 3. Remove 1 mL from test tube T1 and add it to test tube T2. Vortex T2.**
- 4. Remove 1 mL from test tube T2 and add it to test tube T3. Vortex T3.**
- 5. Remove 1 mL from test tube T3 and add it to test tube T4. Vortex T4.**
- 6. Remove 1 mL from test tube T4 and add it to test tube T5. Vortex T5.**
- 7. Remove 1 mL from test tube T5 and add it to test tube T6. Vortex T6.**


Uses of Serial Dilution Method

Serial dilution is one of the most important skills a biology graduate must develop. It has wide applications in various disciplines of biology. It is used for isolating and quantitating the number of microorganisms present in various samples such as **water, food, soil** .

Serial dilution skill is required for bioburden testing

- 1-minimum inhibitory concentration (MIC).**
- 2-most probable number method (MPN).**
- 3-determination of antibody titer.**
- 4-determination of minimum lethal dose (MBC).**

Fill dilution tubes with 90 μL of sterile diluent e.g. SM buffer before you start



Note: change tip and vortex tube each time you serially transfer

