

Bacterial taxonomy

Lab:6

Differential biochemical tests

Indole , Methyl Red , Voges-Proskauer and Citrate utilization tests (IMViC).

Indole Test

The indole test determines a bacterium's ability to produce the enzyme tryptophanase, which breaks down the amino acid tryptophan into indole. It is a key test in identifying Gram-negative bacteria, particularly in the IMViC series (Indole, Methyl red, Voges-Proskauer, Citrate) to differentiate **Enterobacteriaceae**. A positive result is indicated by a **red ring** after adding **Kovac's reagent**, while a negative result remains **yellow**.

Principle and Procedure

- **Principle:** Bacteria possessing the enzyme **tryptophanase** break down tryptophan into **indole**, **ammonia**, and **pyruvic acid**.
- **Reagent:** **Kovac's reagent** (containing p-dimethyl amino benzaldehyde) is added, which reacts with indole to form a red dye (rosin dole).

• Procedure (Tube Method):

1. Inoculate tryptophan or peptone broth with the bacteria.
 2. Incubate for 24–48 hours at 35–37°C.
 3. Add 5 drops of Kovac's reagent.
 4. Observe for a red layer at the top of the broth within minutes.
- **Rapid/Spot Method:** Rub a colony onto filter paper saturated with dimethylamino-cinnamaldehyde (DMACA) reagent. A blue-green color indicates a positive result.

Results

- **Positive (+):** Formation of a red (Kovac's) or blue (spot test) color.
- **Negative (-):** The reagent layer remains yellow or colorless.

Common Bacteria Tested

Indole Positive:

Escherichia coli,

Klebsiella oxytoca, Proteus vulgaris,

Morganella morganii, Vibrio cholerae.

Indole Negative:

Klebsiella pneumoniae


Proteus mirabilis

Salmonella spp

Pseudomonas aeruginosa.

Limitations

- **Media Requirements:** The medium must be rich in tryptophan (e.g., tryptone or peptone broth).
- **Medium Interference:** Do not use indole tests on colonies from MacConkey agar as the medium color can interfere with results.
- **Anaerobes:** For anaerobes, Ehrlich's reagent is preferred over Kovac's reagent



Spot Indole Test

-ve



A

+ve



B

Methyl Red test

The Methyl Red test (MR test) identifies bacteria that perform mixed acid fermentation of glucose, producing stable, acidic end-products, by using the pH indicator methyl red which turns red in acidic conditions (pH below 4.4) and yellow at higher pH. Used with the Voges-Proskauer test (VP) and Indole test (part of IMViC), it helps **differentiate Enterobacteriaceae**, like *E. coli* (positive) from *Enterobacter* (negative). A red color indicates a positive result, while yellow signifies a negative-result.

Principle

- **Bacteria ferment glucose, producing pyruvic acid.**

- **Positive (Mixed Acid Pathway):**

Some bacteria further convert pyruvic acid into significant amounts of stable acids (lactic, acetic, formic, succinic), drastically lowering the pH below 4.4.

- **Negative (Butylene Glycol Pathway):**

Other bacteria convert pyruvic acid to neutral products (like acetoin), resulting in less acidic conditions.

- **Methyl red indicator is added; red means positive (low pH), yellow means negative (higher pH).**

Procedure

- 1. Inoculate MR-VP broth (containing glucose, peptone, buffer) with the test organism.**
- 2. Incubate aerobically at 37°C for 24-48 hours.**
- 3. Add 2-3 drops of methyl red indicator to a portion of the broth.**
- 4. Observe the color change immediately.**


Results

- **Positive: A stable red color.**
- **Negative: A yellow color (or no change from the initial yellow/orange of the broth).**

Common Uses

- **Differentiating bacteria within the Enterobacteriaceae family, such as Escherichia coli (MR positive) from Enterobacter aerogenes (MR negative).**
- **Microbiology identification schemes (IMViC tests).**






Methyl red

Red
+ve MR
(pH approx. 4)

Yellow
-ve MR
(pH approx. 6)



Voges Proskauer test

The **Voges-Proskauer** (VP) test is a biochemical test used to determine if a bacterium produces acetoin (acetyl methyl carbinol) as a neutral end product of glucose fermentation. It is a key test in the IMViC series.

primarily used to differentiate members of the family *Enterobacteriaceae*, such as *Escherichia coli*(VP= -)

from *Klebsiella* and *Enterobacter* species VP=+).

Principle

Bacteria that use the butanediol fermentation pathway metabolize pyruvic acid (a product of glucose metabolism) into acetoin, a neutral intermediate product. The VP test works by detecting this intermediate. When the VP reagents (**alpha-naphthol and potassium hydroxide**) are added, the acetoin in the medium is oxidized to diacetyl in the presence of atmospheric oxygen. Diacetyl then reacts with the guanidine compounds in the peptone broth to form a **pinkish-red polymer**, indicating a positive result. Alpha-naphthol serves as a color intensifier in this reaction.

Procedure

- 1. Inoculation:** A pure culture of the bacterium is aseptically inoculated into a tube containing MR-VP broth (Methyl Red-Voges-Proskauer broth), a buffered peptone-glucose medium.
- 2. Incubation:** The inoculated tube is incubated at 35-37°C for 24 to 48 hours.
- 3. Reagent Addition:** After incubation, the following reagents are added to an aliquot of the culture:
 - 1. 5% alpha-naphthol (VP Reagent A)** is added first.
 - 2. 40% potassium hydroxide (KOH) (VP Reagent B)** is then added.
- 3. The order of reagent addition is critical to prevent false results.**
- 4. Observation:** The tube is gently shaken to expose the medium to atmospheric oxygen and observed for color development for up to **one hour.**



Results and Interpretation

- **VP Positive:** A pink-red or cherry red color develops on the surface of the medium within 15 minutes to an hour, indicating the presence of acetoin.
- **VP Negative:** The medium remains yellow or yellow-brown (or may develop a copper color) after the addition of reagents, indicating the absence of acetoin.

VP Positive Organisms

Specific bacteria are known to be VP positive, including:

- *Klebsiella* species
- *Enterobacter* species
- *Hafnia* species
- *Serratia* species
- *Staphylococcus aureus*
- Most Viridans streptococci (except for some species like *S. mitis*)



Citrate Utilization Test

The citrate utilization test determines if an organism can use sodium citrate as its sole carbon source and ammonium salts as a nitrogen source, typically using Simmons Citrate Agar. A positive result, indicating the enzyme citrase is present, turns the green medium blue (pH > 7.6) due to alkaline byproducts.

common for identifying *Enterobacteriaceae*.

Purpose:

Differentiates gram-negative bacteria based on their ability to metabolize citrate.

Medium: Simmons Citrate Agar

(contains sodium citrate, ammonium salt, and bromothymol blue).

Procedure: A streak-stab technique is used to inoculate a slant, which is incubated at 35-37°C for 24–48 hours.

Interpretation:

- **Positive (+): Growth on the slant and a color change from green to deep blue.**
- **Negative (-): No growth, medium remains green.**

Examples:

***Klebsiella pneumoniae* and *Enterobacter aerogenes* are positive**

while *Escherichia coli* is negative.



**citrate
negative**

**citrate
positive**

