



Biopharmaceutics lab 6 part 2

In vitro dissolution study of per – oral tablet

Procedure (Dissolution study of nitrofurantoin tablets)

- Fill the jars with 1 L of dissolution fluid (artificial gastric juice)
- Put the jars in a thermostatically controlled water bath at 37 (switch water bath on)
- Place 1 tablet of nitrofurantoin in the basket of apparatus, and start the instrument immediately
- Set the speed on 50 r.p.m
- With draw 5 ml each 10 min for 1.5 hr

Notes

- stop apparatus for withdrawal.
- Substitute for the volume withdrawn each time interval using a fresh (artificial gastric juice) previously maintained at 37° C
- After withdrawal and replacement start instrument again
- Analyze samples of nitrofurantoin by reading the absorbance on a spectronic -20 at 369 nm ~370 nm immediately because nitrofurantoin solution discolored by exposure to light
- Use the following straight line equation to obtain the concentration of $y = c + bx$, $y = AB$, $X =$ conc.(mg/100ml), $c = 0.01$, $b = 0.4$




Collection time (min)	Ab	Conc.(mg/ml)	Percent drug release
10			
20			
30			
40			
50			
60			
75			
90			

- To obtain percent drug release:
- $\text{Con. (mg/ml)} \times \text{vol. of dissolution media} = \text{mg amount of drug released into the dissolution media}$
- $\% \text{ of drug released} = \frac{\text{amount of drug release}}{\text{amount of drug in the tablet}} \times 100$
- Plot the following graphs :
- 1. drug release % versus time in (min)
- 2. log % drug released versus time in (min)
- Find dissolution rate K from the slope

$$K = \text{slope} \times 2.303$$

- Note : amount of drug is equal to 100 mg in each tab





WORKING EXAMPLE USING EXCEL SHEET

File Home Insert Page Layout Formulas Data Review View

Themes Colors Fonts Effects Themes Margins Orientation Size Print Area Breaks Background Print Titles Width: Automatic Height: Automatic Scale: 100% Sheet Options Gridlines View Print View Print Bring Forward Send Backward Selection Pane Align Group Rotate Arrange

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R
1																		
2																		
3																		
4																		
5																		
6																		
7																		
8																		
9																		
10																		
11																		
12																		
13																		
14																		
15																		
16																		
17																		
18																		
19																		
20																		
21																		
22																		
23																		
24																		
25																		

Switch the sheet direction so that the first column is on right side

Ready 11:10 AM 12/3/2018

Type here to search



absorbance	conc u/ml	amount (m	% drug released
0	0	0	0
0.165	$=B3/0.0195$		
0.223			
0.327			
0.283			

in vitro dissolution test was done t
Drug A (Dose 10 mg)
dissolution medium volume was 50
I.D.W.
libration curve is
0.0195 X

First Reading is always zero

calculate Conc $\mu\text{g/ml}$ with the help of $Y = bx + c$ convert this equation in $x = Y - c / m$

Always starting calculation with =

Note units of x are given in u/ml in this example



Microsoft Excel ribbon: File, Home, Insert, Page Layout, Formulas, Data, Review, View

Font: Calibri, 11, Bold, Italic, Underline, Color, Background Color, Text Color, Text Direction, Paragraph Spacing, Bullets, Numbering, Indentation, Merge & Center

Number: General, Currency, Percentage, Decimals, Thousands Separator, Comma Separator, Text to Columns, Text to Rows

Styles: Conditional Formatting, Format as Table, Cell Styles

Cells: Insert, Delete, Format

Editing: AutoSum, Fill, Clear, Sort & Filter, Find & Select

1	Time(minut	absorbance	conc u/ml)	amount (m	% drug released
2	0	0	0	0	0
3	60	0.165	8.461538		
4	120	0.223	11.4359		
5	180	0.327	16.76923		
6	240	0.283	14.51282		

Formula bar: C3 = (B3/0.0195)

Cell C3: 8.461538

Cell C4: 11.4359

Cell C5: 16.76923

Cell C6: 14.51282

Statistical summary: Average: 12.79487179, Count: 4, Sum: 51.17948718

System tray: 12:42 PM, 12/3/2018

After completing calculation press ENTER Then DRAG

in vitro dissolution test was done to Drug A (Dose 10 mg)
Dissolution medium volume was 500 ml D.W.
calibration curve is $y=0.0195X$

Time (min)	absorbance	conc u/ml	amount (mg)	% drug released
0	0	0	0	0
60	0.165	8.461538	=(C3*500)/1000	
120	0.223	11.4359		
180	0.327	16.76923		
240	0.283	14.51282		

amount in mg = conc u/ml * volume of media / 1000

in vitro dissolution test was done to Drug A (Dose 10 mg)
Dissolution medium volume was 500 ml D.W.
calibration curve is y=0.0195 X



Microsoft Excel interface showing a spreadsheet with columns for Time, absorbance, concentration, amount, and % drug released. A formula $= (D3/10) * 100$ is used to calculate the percentage of drug released. A callout box explains the formula: Percent drug released = (amount released/amount in one tablet*100). A text box provides experimental details: in vitro dissolution test was done to Drug A (Dose 10 mg), Dissolution medium volume was 500 ml D.W., calibration curve is $y=0.0195X$.

Time (min)	absorbance	conc (u/ml)	amount (mg)	% drug released
0	0	0	0	0
60	0.165	8.461538	4.230769	$= (D3/10) * 100$
120	0.223	11.4359	5.717949	
180	0.327	16.76923	8.384615	
240	0.283	14.51282	7.25641	

Percent drug released = (amount released/amount in one tablet*100)

in vitro dissolution test was done to Drug A (Dose 10 mg)
Dissolution medium volume was 500 ml D.W.
calibration curve is $y=0.0195X$



example lab 5 - Microsoft Excel

File Home Insert Page Layout Formulas Data Review View

PivotTable Table Picture Clip Art Shapes SmartArt Screenshot Column Line Pie Bar Area Scatter Other Charts Line Column Win/Loss Slicer Hyperlink Text Box Header & Footer WordArt Signature Line Object Equation Symbol

Tables Illustrations Charts Sparklines Filter Links Text Symbols

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R
1	Time(minut	absorbance	conc u/ml)	amount (m	% drug released													
2	0	0	0	0	0													
3	60	0.165	8.461538	4.230769	42.30769													
4	120	0.223	11.4359	5.717949	57.17949													
5	180	0.327	16.76923	8.384615	83.84615													
6	240	0.283	14.51282	7.25641	72.5641													
7																		
8	Time(minut	% drug released																
9	0	0																
10	60	42.30769																
11	120	57.17949																
12	180	83.84615																
13	240	72.5641																

in vitro dissolution test was done to Drug A (Dose 10 mg)
Dissolution medium volume was 500 ml D.W.
calibration curve is $y=0.0195X$

Scatter with Smooth Lines and Markers
Compare pairs of values.
Use it when there are a few data points in x-axis order and the data represents a function.

All Chart Types...

1. Select (time and % drug release columns)
2. Insert
3. Scatter with smooth lines
4. enter

Select destination and press ENTER or choose Paste

Average: 85.58974359 Count: 12 Sum: 855.8974359 100%

Type here to search

12:51 PM 12/3/2018



example lab 5 - Microsoft Excel

ChartTools: Design, Layout, Format

Chart Name: Chart 2

Vertical (Value) Axis Major Gridlines

Primary Horizontal Axis Title

Primary Vertical Axis Title

To add titles

Time (min)	absorbance	conc u/ml	amount (mg)	% drug released
0	0	0	0	0
60	0.165	8.461538	4.230769	42.30769
120	0.223	11.4359	5.717949	57.17949
180	0.327	16.76923	8.384615	83.84615
240	0.283	14.51282	7.25641	72.5641

Time (min) % drug released

0	0
60	42.30769
120	57.17949
180	83.84615
240	72.5641

in vitro dissolution test was done to Drug A (Dose 10 mg)
Dissolution medium volume was 500 ml D.W.
calibration curve is $y=0.0195X$

Ready


12:53 PM 12/3/2018



File Home Insert Page Layout Formulas Data Review View

Clipboard Font Alignment Number Styles Cells Editing

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R
1	Time(minut	absorbance	conc u/ml)	amount (m	% drug released													
2	0	0	0	0	0													
3	60	0.165	8.461538	4.230769	42.30769													
4	120	0.223	11.4359	5.717949	57.17949													
5	180	0.327	16.76923	8.384615	83.84615													
6	240	0.283	14.51282	7.25641	72.5641													
7																		
8	Time(minut	% drug released																
9	0	0																
10	60	42.30769																
11	120	57.17949																
12	180	83.84615																
13	240	72.5641																
14																		
15																		
16																		
17																		
18																		
19																		
20																		
21																		
22																		
23																		
24																		
25																		



in vitro dissolution test was done to Drug A (Dose 10 mg)
Dissolution medium volume was 500 ml D.W.
calibration curve is $y=0.0195X$

143,9 243,9 343,9

Enter 100% 12:58 PM 12/3/2018



Thank you