Mustansiriyah University
College Of Science
Department Of Chemistry

ORGANIC IDENTIFICATION LABORATORY

Asmaa A. Al-Mutash

Determination of functional group and general test

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The determination of the functional group depends largely on the correct determination of elements once the functional group present in the organic compound is known, one can find out the name of the probable compound with the help of m.p and p.b. next procedure is to study some specific reactions of the compound and then to confirm its name by preparing suitable derivatives. It may be noted that elements C, H, and O are not tested. This means that functional groups containing C, H, and O are determined if no other element exists.

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- A. Oxygen organic Compounds.
- B. Nitrogen Compounds.
- C. Aromatic Hydrocarbon Compounds.
- D. Alkyl Halide and Aryl halide Compounds.

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
a-litmus test: -

Shake a small amount of the given substance in water and immerse a blue litmus paper (L.p.) in it. If the L.p. turns **red** means presence of Carboxylic acid (phenols being acidic in nature also give this test).

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**b- Sodium bicarbonate test: -**



Add a pinch of the given organic substance if solid (2 drops of liquid) to 5 ml of saturated solution (aqueous) of NaHCO_3 and shake gently. The effectiveness or slow up and down movement of particles indicates the presence of a carboxylic group.

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
Phenols < carbonic acid < carboxylic acid

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a-litmus test: -

Blue litmus paper change to **red** but no effervescence with NaHCO₃.


b- FeCl₃ test: -

Dissolve the organic substance in water in clean test tube and add drop of aq. Solution of FeCl₃.

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






To 1 ml of organic substance in H_2O , add a few drops of ceric ammonium nitrate change the yellow color of indicator to red that is means alcohol is present.

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




b- Chromic acid test. (Jones oxidation).

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
Primary alcohol $\longrightarrow$ Aldehyde.

Secondary alcohol $\longrightarrow$ Ketone.

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







Procedure:

In clean test tube dissolve 2 drops of alcohol in acetone and add 2 drops from chromic acid reagent change the color to green or blue green means presence (1°, 2° alcohols).

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

Procedure: -

Add 5 ml from Lucas test to 1 ml of organic substant in clean and dry test tube, shake well and notice: -

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








a- 2,4-dinitro phenyl hydrazine (brady Test): -

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
Procedure: -

Add (1 mL) 2,4-dinitro phenyl hydrazine (2,4-DNPH) reagent to organic substance in water or Methanol (**safety note:** the reagent is highly toxic) and mix the test tube by agitating. we see colored ppt appear.

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
The Tollen's reagent is a mild oxidizing agent that can oxidize aldehydes, but not alcohols or other carbonyl compounds. A positive test result is the formation of the element silver which precipitates out as a "silver mirror" on the test tube, or as a black colloidal precipitate.

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






(left to right) 2-pentanone (ketone)(-ve), isobutyraldehyde, benzaldehyde (+ve)

General equation: -**C – Iodoform Test: - KI /I₂**

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
**Procedure: -**

add KI /I₂ (3drops) to organic substance in 10% NaOH, formation of ppt indicate presence of $-\text{CH}_3$ beside carbonyl group.

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* a) Iodoform results (left to right): 1-propanol, 2-propanol, acetone, 1-octanol (b) Negative result, c) Positive result, d) Negative result when the sample is insoluble in the reagent.

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
Note:

*A positive result is a cloudy yellow solution, or a yellow precipitate. A negative result is a clear, yellow, or orange solution with no precipitate.


Note:

a false positive result may occur if the test tube was cleaned with acetone before use, and residual acetone remained in the tube.

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


Procedure: -

Dissolve 0.5 gm of organic substance in 2 ml methanol. To it add 4-5 drops dil. NaOH solution and 2-3 drops of ph. indicator. A pink color is produced. Now put the test tube in boiling water bath for 5-10 mint. If the pink color either fades or disappears, this shows presence of ester group.

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
Procedure: -

Take 2 drops of Unknown or few grams of organic substance in a test. tube contain methanol, add 1ml of 5% methanolic solution of hydroxyl amine hydrochloride. Now, add a few drops of (2N) Methanolic KOH until the solution is just alkaline, Boil, Cool and acidify with dil. HCl until the solution became neutral (زيادة الحامض غير مفيدة), add FeCl_3 . If deep red or violet color produced, it indicates the presence of ester group.


a) Preliminary procedure to be sure the test will not give a false positive, b) Heating the solution, c) Negative result (2-pentanone), d) Positive result (ethyl acetate).

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procedure: -

To 2ml aq. Solution of organic substance add (2-3 drops) of Molisch's reagent (alcoholic ' α - naphthol) and shake well. Now add carefully (1 ml) of conc. H_2SO_4 from the side of the test tube. Allow to stand. A reddish violet ring of the junction of two liquids shows the presence of carbohydrate.

**b – Blackening with H_2SO_4**

Heat 0.2 gm of organic substance with 1ml of conc. H_2SO_4 . An immediate charring and blackening of the solution show presence of carbohydrate.

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